cYmedics™

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
The screenshot displays the 'patients database' interface of the 'cymedics yaoteq' software, version 2012. The interface is dark-themed with orange and blue accents. At the top, there are navigation buttons for 'Support', 'Settings', and 'Exit'. Below these, there are dropdown menus for 'acupuncture mode' (set to 'Acupuncture') and 'instrument', along with 'Rescan' and 'Reset' buttons. A left sidebar contains several menu items: 'Patients Database' (selected), 'Manage Patient's Clinical Information', 'Manage Project Point Data Sets', 'Acquire Data From Patient', 'Analyze Patient's Data', and 'Operate Patient's Energy Correction'. The main area shows a search bar with the text 'patients database' and a 'New' button. Below the search bar is a grid of letters A-Z. A table below the grid lists patient information:

Surname	First Name	Birthdate	Sex
Platz	Sabine	03/04/1966	Female

At the bottom of the screen, there is a status bar showing '[No Patient]', 'status Ready.', and environmental data: 'air 23 °C / 70 %RH' and '27.02.2012'.

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
The screenshot shows the 'acupuncture data acquisition' screen of the 'cymedics yaoteq' software, version 2012. The interface is dark-themed with orange and blue accents. At the top, there are buttons for 'Support', 'Settings', and 'Exit'. Below that, a dropdown menu shows 'acupuncture mode' set to 'Acupuncture', and another dropdown shows 'instrument' with 'Rescan' and 'Reset' buttons. On the left, a sidebar contains navigation options: 'Patients Database', 'Manage Patient's Clinical Information', 'Manage Project', 'Point Data Sets', 'Acquire Data From Patient', 'Analyze Patient's Data', and 'Operate Patient's Energy Correction'. The main area is titled 'acupuncture data acquisition' and features a 'point set' dropdown set to 'Standard' and a 'Single' checkbox. Below this is a table with two columns: 'Point' and 'Potential'. The table lists several points with a potential of 0. To the right of the table is an image of a human arm with a white line tracing the meridian path and numbered points 1 through 5. At the bottom of the main area, there is a 'Start' button and the text 'Press Start to begin the acquisition.' The bottom status bar shows 'Platz, Sabine' as the user, 'Ready' as the status, and environmental data: '23 °C / 70 %RH' and the date '27.02.2017'.

Point	Potential
Yg Lu 5 Chize	0
Yg Lu 8 Jingqu	0
Yg Lu 9 Taiyuan	0
Yg Lu 10 Yuji	0
Yg Lu 11 Shaoshang	0
Yg LI 1 Shangyang	0
Yg LI 2 Erjian	0

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The screenshot shows the 'acupuncture data acquisition' screen of the 'cymedics yaoteq version 2012' software. The interface includes a top navigation bar with 'Support', 'Settings', and 'Exit' buttons. Below this, there are controls for 'acupuncture mode' (set to 'Acupuncture') and 'instrument' (with 'Rescan' and 'Reset' buttons). A left sidebar contains menu items: 'Patients Database', 'Manage Patient's Clinical Information', 'Manage Project', 'Point Data Sets', 'Acquire Data From Patient' (highlighted), 'Analyze Patient's Data', and 'Operate Patient's Energy Correction'. The main area is titled 'acupuncture data acquisition' and features a 'point set' dropdown (set to 'Standard') and a 'Single' checkbox. Below these is a table of points with their potentials. To the right, there are 'image type' options: 'Point' (checked), 'Meridian', and 'Pathology'. A 3D anatomical model of a human leg and foot shows the meridian path with numbered points. At the bottom, there is a 'Start' button and the instruction 'Press Start to begin the acquisition.' The status bar at the very bottom shows patient information: 'Platz, Sabine', 'status: Ready', '23 °C / 70 %RH', and '27.02.2017'.

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The screenshot displays the 'cymedics yaoteq' software interface, version 2012. The main window is titled 'acupuncture data analysis'. On the left, there is a sidebar with a 'Patients Database' icon and several menu items: 'Manage Patient's Clinical Information', 'Manage Project', 'Point Data Sets', 'Acquire Data From Patient', 'Analyze Patient's Data', and 'Operate Patient's Energy Correction'. The main area shows a calendar for February 2013 with the 27th highlighted. Below the calendar is a table with columns 'Point' and 'Potential'. The table is currently empty. At the bottom of the interface, there is a status bar showing 'Platz, Sabine' with a 'Ready' status, and environmental data: '23 °C / 70 %RH' and the date '27.02.2013'.

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Bibliografia

- **Napadow V, Makris N, Liu J, Kettner N, Kwong K, Hui K**, (2005) Effects of electroacupuncture versus Manual Acupuncture on the Human Brain as Measured by fMRI.
- **Pischinger A** (2004) Das System der Grundregulation – Grundlagen einer ganzheitsbiologischen Medizin. 10. Durchges. Aufl., bearbeitet von Hartmut Heine, Otto Bergsmann und Felix Perger. Haug, Stuttgart
- **Starwynn D**, (2006) Is Microcurrent Therapy Effective for Treatment of Macular Degeneration? Acupuncture Today, May 2006, Vol. 07, Issue 05
- **Matsumoto K, Birch S**, (1988) Hara diagnosis: reflections on the sea. Paradigm, Brookline
- **Morton M A, Dlouhy C** (1989) Energy fields in medicine. John E. Fetzer Foundation, Kalamazoo, MI
- **Oschmann JL, Oschman NH** (1995) Physiological and emotional effects of acupuncture needle insertion. Proceedings of the Second Symposium of the Society for Acupuncture Research, SAR, Boston
- **Reichmanis M, Marino A A, Becker R O** (1975) Electrical correlates of acupuncture points. IEEE Transactions on Biomedical Engineering, 22 (November): 533-535
- **Becker R O** (1990) The machine brain and properties of the mind. Subtle Energies 113:79-97
- **Becker R O** (1991) Evidence for a primitive DC electrical analog system controlling brain function. Subtle Energies 2(1):71-88
- **Fröhlich H** (1968) Bose condensation of strongly excited longitudinal electric modes. Physics Letters 26A:402-403
- **Fröhlich H, ed** (1988) Biological coherence and response to external stimuli. Springer Verlag, Berlin
- **Heidemann SR** (1993) A new twist on integrins and the cytoskeleton. Science 260:1080-1081
- **Ho M-W, Knight DP** 1998 The acupuncture system and the liquid crystalline collagen fibers of the connective tissues. American Journal of Chinese Medicine 26(3-4):1-13
- **Ingber DE** (1993a) The riddle of morphogenesis: a question of solution chemistry or molecular cell engineering? Cell 75:1249-1252
- **Ingber DE** (1993b) Cellular tensegrity: defining new rules of biological design that govern the cytoskeleton. Journal of Cell Science 104:613-627
- **Ingber DE** (1998) The architecture of life. Scientific American 278(1):48-57
- **Oschmann JL** (1994) A biophysical basis for acupuncture. Proceedings of the First Symposium of the Committee for Acupuncture Research, CAR, Boston
- **Oschmann JL** (1996) The nuclear, cytoskeletal, and extracellular matrices: a continuous communication network. Poster presentation for “the cytoskeleton: mechanical, physical and biological interactions”, a workshop sponsored by the Center for Advanced Studies in the Space Life Science at the marine Biological Laboratory, Woods Hole, Massachusetts, November 15-17, 1996
- **Oschman JL, Oschman NH** (1994) Book review and commentary: Biological coherence and response to external stimuli, edited by Herbert Fröhlich, published by Springer Verlag, Berlin, 1988. NORA Press, Dover, NH
- **Pienta KJ, Coffey DS** (1991) Cellular harmonic information transfer through a tissue tensegrity-matrix system. Medical Hypotheses 34:88-95
- **Pugh A** (1976) An introduction to tensegrity. University of California Press, Berkeley
- **Wang JY, Butler JP, Ingber DE** (1993) Mechanotransduction across the cell surface and through the cytoskeleton. Science 260:1124-1127